

# WORKFLOWS FOR PFAS TESTING



## PerkinElmer PFAS QSight® LC/MS Quantitative Workflow


Adding new capabilities to your analytical laboratory in an efficient way can be challenging – especially without predetermined methods. So to accelerate this process, we've developed workflows for EPA PFAS methods, including Method 533, Method 537.1, and Draft Method 1633. Once these methods are established in your busy production environment, maximizing sample throughput is key. Our innovative StayClean technology enables you to run more samples per week due to a robust tolerance and minimal cleaning requirements for dirty matrices such as sludge, tissue, food, and soil.



### INCREASED ROBUSTNESS AND UPTIME

Revolutionary StayClean/HSID technologies allow you to analyze dirty, complex matrices. Here's how:

- Continuous flow of hot gas acts as a constant cleaning agent
- Entirely flow dominant, it requires no optimization or voltage adjustment




## PFAS Workflow at a Glance



# PFAS Workflow = Fastest Runtime, Maximum Uptime

## EPA DRAFT METHOD 1633



## QSight EPA Draft Method 1633 Performance

Analyte	Correlation Coefficient (R <sup>2</sup> )	LOD (ng/L)	LOQ (ng/L)
PFBA	0.999941	0.2	0.6
PFMPA	0.999993	0.09	0.3
PFPeA	0.999952	0.04	0.1
3:3FTCA	0.999989	1	5
PFBS	0.999979	0.2	0.5
PFMBA	0.999991	0.06	0.2
PFEESA	0.999977	0.03	0.1
NFDHA	0.999880	0.08	0.3
4:2FTS	0.999921	0.09	0.3
PFHxA	0.999931	0.04	1
PFPeS	0.999970	0.03	0.1
HFPO-DA	0.999962	0.4	1
PFHxS	0.999880	0.03	0.1
PFHpA	0.999974	0.06	0.2
ADONA	0.999585	0.01	0.1
5:3FTCA	0.999388	0.3	1
6:2FTS	0.998310	0.1	0.4
PFOA	0.999987	0.1	0.5
PFHpS	0.999868	0.02	0.1
PFOS	0.999950	0.03	0.1

Analyte	Correlation Coefficient (R <sup>2</sup> )	LOD (ng/L)	LOQ (ng/L)
PFNA	0.999721	0.2	0.6
7:3FTCA	0.999784	0.3	1
9Cl-PF3ONS	0.999850	0.5	2
NMeFOSAA	0.999326	0.02	0.1
PFNS	0.999903	0.2	0.5
PFDA	0.999913	0.5	2
8:2FTS	0.999629	0.1	0.4
NEtFOSAA	0.999965	0.1	0.4
PFDS	0.999901	0.04	0.1
PFUnA	0.999970	0.5	2
11-PF3OUds	0.999747	0.01	0.1
PFOSA	0.999730	0.01	0.05
PFDoA	0.999919	0.3	1
PFTDA	0.999724	0.2	0.8
PFDoS	0.999899	0.5	2
PFTeDA	0.999661	0.2	0.6
NMeFOSA	0.999976	0.1	0.4
NMeiFOSE	0.999986	0.9	3
NEtFOSE	0.999978	0.3	0.9
NEtFOSA	0.999963	0.07	0.2

EPA Method 537.1 Compounds

ORDER YOUR FREE PFAS POSTER



PerkinElmer, U.S. LLC.  
710 Bridgeport Ave.  
Shelton, CT 06484-4794 USA  
(+1) 855-726-9377  
[www.perkinelmer.com](http://www.perkinelmer.com)



For a complete listing of our global offices, visit [www.perkinelmer.com/ContactUs](http://www.perkinelmer.com/ContactUs)

Copyright ©2023, PerkinElmer, U.S. LLC. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.